

# AST-1000

Luxury Style  
Ophthalmic Unit Table

## • Overall

- Unit Size 2,000 X 1,150 X 1,730 mm
- Input Voltage 110~220 V 50/60 Hz

## • Chair

- Electric Up-Down Movement (510~660 mm)
- Electric Front-Rear Movement (100 mm)

## • Electric Arm



- Electric Arm Lifting Movement  
Up-Down 150 mm  
Test-Rest 400 mm
- Near Chart LED Lamp 3 W

## • Sliding Table



- Load 2 Devices
- Dimension 1,250 X 460 mm
- Electric Up-Down Movement (770~950 mm)
- Electric Sliding Movement (560 mm)

## • Drawer



- Inner Size  
470 X 360 X 50 mm
- 2 Drawers



# AST-1000

MODEL	COMPONENT	WEIGHT	CBM	PACKAGE
AST-1000	CABINET	300.00	1.96	WOOD
	CHAIR	60.00	0.37	PAPER
HRT-1	CABINET+TABLE	200.00	0.92	WOOD
	CHAIR	70.00	0.40	PAPER
HRT-1A	CABINET+TABLE	200.00	0.92	WOOD
	CHAIR	70.00	0.40	PAPER
RT-9000	TABLE	35.00	0.19	PAPER
	CABINET	160.00	0.63	PLYWOOD
	CHAIR	45.00	0.40	PAPER
C11	TABLE	30.00	0.17	PAPER
	CABINET	150.00	0.63	PLYWOOD
	CHAIR	45.00	0.40	PAPER
C10	TABLE	35.00	0.19	PAPER
	CABINET	160.00	0.63	PLYWOOD
	CHAIR	45.00	0.40	PAPER
C6	TABLE	30.00	0.17	PAPER
	CABINET	150.00	0.63	PLYWOOD
	CHAIR	38.00	0.34	PAPER
TOTAL		1,823.00	9.75	