

HRT-1A

Luxury Style
Ophthalmic Unit Table

- Overall

- Unit Size 1,470 X 1,082 X 2,025 mm
- Input Voltage 110~220 V 50/60 Hz

- Chair

- Electric Up-Down Movement (550~650 mm)
- Electric Front-Rear Movement (100 mm)

- LED Lamp 3 W

- Swing Arm

- Swing Arm Lifting Movement
Up-Down 40 Degree
Right-Left 180 Degree

- Sliding Table

- Load 2 Devices
- Dimension 980 X 386 mm
- Rotation Angle 90 Degree
- Electric Up-Down Movement (770~870 mm)
- Electric Sliding Movement (340 mm)

- Drawer

- Inner Size
510 X 360 X 55 mm



HRT-1A

MODEL	COMPONENT	WEIGHT	CBM	PACKAGE
AST-1000	CABINET	300.00	1.96	WOOD
	CHAIR	60.00	0.37	PAPER
HRT-1	CABINET+TABLE	200.00	0.92	WOOD
	CHAIR	70.00	0.40	PAPER
HRT-1A	CABINET+TABLE	200.00	0.92	WOOD
	CHAIR	70.00	0.40	PAPER
RT-9000	TABLE	35.00	0.19	PAPER
	CABINET	160.00	0.63	PLYWOOD
	CHAIR	45.00	0.40	PAPER
C11	TABLE	30.00	0.17	PAPER
	CABINET	150.00	0.63	PLYWOOD
	CHAIR	45.00	0.40	PAPER
C10	TABLE	35.00	0.19	PAPER
	CABINET	160.00	0.63	PLYWOOD
	CHAIR	45.00	0.40	PAPER
C6	TABLE	30.00	0.17	PAPER
	CABINET	150.00	0.63	PLYWOOD
	CHAIR	38.00	0.34	PAPER
TOTAL		1,823.00	9.75	