

HRT-3

Luxury Style
Ophthalmic Unit Table

• Overall

- Unit Size 1,800 X 1,500 X 1,900 mm
- Input Voltage 110~220 V 50/60 Hz

• Chair



- 180 Degree Backrest Fold
- Electric Up-Down Movement (550~650 mm)
- Electric Front-Rear Movement (100 mm)



• Main Lamp 20 W

• Electric Arm



- Electric Arm Lifting Movement
Up-Down 150 mm
Test-Rest 400 mm
- Near Chart Lamp 3 W

• Sliding Table



- Load 2 Devices
- Dimension 965 X 385 mm
- Rotation Angle 90 Degree
- Electric Up-Down Movement (810~900 mm)
- Electric Sliding Movement (340 mm)

• Drawer



- Inner Size
500 X 360 X 80 mm
500 X 360 X 120 mm

HRT-3

MODEL	COMPONENT	WEIGHT	CBM	PACKAGE
AST-1000	CABINET	300.00	1.96	WOOD
	CHAIR	60.00	0.37	PAPER
HRT-1	CABINET+TABLE	200.00	0.92	WOOD
	CHAIR	70.00	0.40	PAPER
HRT-1A	CABINET+TABLE	200.00	0.92	WOOD
	CHAIR	70.00	0.40	PAPER
RT-9000	TABLE	35.00	0.19	PAPER
	CABINET	160.00	0.63	PLYWOOD
	CHAIR	45.00	0.40	PAPER
C11	TABLE	30.00	0.17	PAPER
	CABINET	150.00	0.63	PLYWOOD
	CHAIR	45.00	0.40	PAPER
C10	TABLE	35.00	0.19	PAPER
	CABINET	160.00	0.63	PLYWOOD
	CHAIR	45.00	0.40	PAPER
C6	TABLE	30.00	0.17	PAPER
	CABINET	150.00	0.63	PLYWOOD
	CHAIR	38.00	0.34	PAPER
TOTAL		1,823.00	9.75	